



For the integration of the GeniX system and its accessories such as the water to air chiller system, Xenocs proposes a 19" cabinet.

This cabinet is compact and perfectly suited for home laboratory environment offering accommodation and protection of your electronics, to EMC and IP levels (up to 20).


Advantages to the user

- Door and rear panel combination freely selectable
- All uprights with 19" cut-outs and EIA grid holes
- Inset cover
- Colour combination RAL 7021 / RAL 7035
- Static load-carrying capacity 200 kg

Subject to technical changes without notice


DELIVERY QUANTITY

ITEM	DESCRIPTION	QUANTITY
■ 1	19" frame with 19" cut-outs and EIA grid holes, Al die-cast frame, Al uprights, RAL 7021	1
■ 2	Side panels, St, 1 mm, RAL 7035	2
■ 3	Flush cover, St, 1 mm, RAL 7021	1
■ 4	Foot (up to 16 U made of rubber, from 20 U made of Al)	4



TECHNICAL DATA

■ Usable height	U	12
■ Height (H)	mm	589
■ Width (W)	mm	553
■ Depth (D)	mm	600
■ C	mm	507
■ Weight	Kg	20



DMC-060705-GeniX accessories-19" RACK-01